








AlloDerm SELECT™ Regenerative Tissue Matrix

Product Description	Piece Size (cm)	Coverage (cm ²)	Thickness	Thickness Range (mm)	SKU	<input type="checkbox"/>	QTY	
AlloDerm SELECT™ Regenerative Tissue Matrix  Contour Perforated	X-Large 11.8 x 23.7	200	Medium	1.6 ± 0.4	CXL1518P	<input type="checkbox"/>		
			Thick	2.4 ± 0.4	CXL1519P	<input type="checkbox"/>		
	Large 10.7 x 21.5	164	Thin	1.0 ± 0.2	CL1516P	<input type="checkbox"/>		
			Medium	1.6 ± 0.4	CL1518P	<input type="checkbox"/>		
	Thick	2.4 ± 0.4	132	Thin	1.0 ± 0.2	CM1516P	<input type="checkbox"/>	
				Medium	1.6 ± 0.4	CM1518P	<input type="checkbox"/>	
	Thick	2.4 ± 0.4	77	Thin	1.0 ± 0.2	CS1516P	<input type="checkbox"/>	
				Medium	1.6 ± 0.4	CS1518P	<input type="checkbox"/>	
	Thick	2.4 ± 0.4	77	Thick	2.4 ± 0.4	CS1519P	<input type="checkbox"/>	
				Thick	2.4 ± 0.4	CS1519P	<input type="checkbox"/>	
	AlloDerm SELECT™ Regenerative Tissue Matrix  Contour	X-Large 11.8 x 23.7	200	Medium	1.6 ± 0.4	CXL1518	<input type="checkbox"/>	
				Thick	2.4 ± 0.4	CXL1519	<input type="checkbox"/>	
Large 10.7 x 21.5		164	Thin	1.0 ± 0.2	CL1516	<input type="checkbox"/>		
			Medium	1.6 ± 0.4	CL1518	<input type="checkbox"/>		
Thick		2.4 ± 0.4	132	Thin	1.0 ± 0.2	CM1516	<input type="checkbox"/>	
				Medium	1.6 ± 0.4	CM1518	<input type="checkbox"/>	
Thick		2.4 ± 0.4	77	Thin	1.0 ± 0.2	CS1516	<input type="checkbox"/>	
				Medium	1.6 ± 0.4	CS1518	<input type="checkbox"/>	
Thick		2.4 ± 0.4	77	Thick	2.4 ± 0.4	CS1519	<input type="checkbox"/>	
				Thick	2.4 ± 0.4	CS1519	<input type="checkbox"/>	
AlloDerm SELECT RESTORE™ Regenerative Tissue Matrix  Perforated		Large 20.1 x 23.6	327	Medium	1.6 ± 0.4	RL1518P	<input type="checkbox"/>	
				Thick	2.4 ± 0.4	RL1519P	<input type="checkbox"/>	
	X-Thick			3.4 ± 0.6	RL1522P	<input type="checkbox"/>		
AlloDerm SELECT RESTORE™ Regenerative Tissue Matrix  	Large 20.1 x 23.6	327	Medium	1.6 ± 0.4	RL1518	<input type="checkbox"/>		
			Thick	2.4 ± 0.4	RL1519	<input type="checkbox"/>		
			X-Thick	3.4 ± 0.6	RL1522	<input type="checkbox"/>		





AlloDerm SELECT™ Regenerative Tissue Matrix

Product Description	Piece Size (cm)	Coverage (cm ²)	Thickness	Thickness Range (mm)	SKU	<input type="checkbox"/>	QTY
AlloDerm SELECT™ Regenerative Tissue Matrix  Perforated Rectangle	16 x 20	320	Medium	1.6 ± 0.4	1518320P	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	1519320P	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	1522320P	<input type="checkbox"/>	
AlloDerm SELECT™ Regenerative Tissue Matrix  Rectangle	16 x 20	320	Medium	1.6 ± 0.4	1518320	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	1519320	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	1522320	<input type="checkbox"/>	
	12 x 20	240	Medium	1.6 ± 0.4	15181220	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	15191220	<input type="checkbox"/>	
	12 x 16	192	Medium	1.6 ± 0.4	15181216	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	15191216	<input type="checkbox"/>	
	12 x 12	144	Medium	1.6 ± 0.4	1418144	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	1419144	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	1422144	<input type="checkbox"/>	
	10 x 20	200	Medium	1.6 ± 0.4	15181020	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	15191020	<input type="checkbox"/>	
	10 x 16	160	Medium	1.6 ± 0.4	15181016	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	15191016	<input type="checkbox"/>	
	8 x 20	160	Medium	1.6 ± 0.4	1518160	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	1519160	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	1522160	<input type="checkbox"/>	
	8 x 16	128	Thin	1.0 ± 0.2	1516128	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	1518128	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	1519128	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	1522128	<input type="checkbox"/>	
	6 x 18	108	Medium	1.6 ± 0.4	1518108	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	1519108	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	1522108	<input type="checkbox"/>	
6 x 16	96	Thin	1.0 ± 0.2	1516616	<input type="checkbox"/>		
		Medium	1.6 ± 0.4	1518616	<input type="checkbox"/>		
		Thick	2.4 ± 0.4	1519616	<input type="checkbox"/>		
		X-Thick	3.4 ± 0.6	1522616	<input type="checkbox"/>		

AlloDerm SELECT™ Regenerative Tissue Matrix

Product Description	Piece Size (cm)	Coverage (cm ²)	Thickness	Thickness Range (mm)	SKU	<input type="checkbox"/>	QTY
AlloDerm SELECT™ Regenerative Tissue Matrix  Rectangle	8 x 12	96	Medium	1.6 ± 0.4	151896	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	151996	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	152296	<input type="checkbox"/>	
	6 x 12	72	Thin	1.0 ± 0.2	151672	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	151872	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	151972	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	152272	<input type="checkbox"/>	
	4 x 16	64	Thin	1.0 ± 0.2	151664	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	151864	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	151964	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	152264	<input type="checkbox"/>	
	5 x 10	50	Thin	1.0 ± 0.2	141650	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	141850	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	141950	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	142250	<input type="checkbox"/>	
	4 x 12	48	X-Thin	0.55 ± 0.25	141448	<input type="checkbox"/>	
			Thin	1.0 ± 0.2	151648	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	151848	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	151948	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	152248	<input type="checkbox"/>	
	4 x 7	28	X-Thin	0.55 ± 0.25	141428	<input type="checkbox"/>	
			Thin	1.0 ± 0.2	141628	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	141828	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	141928	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	142228	<input type="checkbox"/>	
	2 x 12	24	X-Thick	3.4 ± 0.6	142224	<input type="checkbox"/>	
	3 x 7	21	Thin	1.0 ± 0.2	141621	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	141821	<input type="checkbox"/>	
Thick			2.4 ± 0.4	141921	<input type="checkbox"/>		
X-Thick			3.4 ± 0.6	142221	<input type="checkbox"/>		
2 x 4	8	X-Thin	0.55 ± 0.25	141408	<input type="checkbox"/>		
		Thin	1.0 ± 0.2	141608	<input type="checkbox"/>		
		Medium	1.6 ± 0.4	141808	<input type="checkbox"/>		
		Thick	2.4 ± 0.4	141908	<input type="checkbox"/>		
		X-Thick	3.4 ± 0.6	142208	<input type="checkbox"/>		

AlloDerm SELECT™ Regenerative Tissue Matrix

Product Description	Piece Size (cm)	Coverage (cm ²)	Thickness	Thickness Range (mm)	SKU	<input type="checkbox"/>	QTY
AlloDerm SELECT™ Regenerative Tissue Matrix  Rectangle	1 x 4	4	X-Thin	0.55 ± 0.25	141404	<input type="checkbox"/>	
			Thin	1.0 ± 0.2	141604	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	141804	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	141904	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	142204	<input type="checkbox"/>	
	1 x 2	2	X-Thin	0.55 ± 0.25	141402	<input type="checkbox"/>	
			Thin	1.0 ± 0.2	141602	<input type="checkbox"/>	
			Medium	1.6 ± 0.4	141802	<input type="checkbox"/>	
			Thick	2.4 ± 0.4	141902	<input type="checkbox"/>	
			X-Thick	3.4 ± 0.6	142202	<input type="checkbox"/>	
	1 x 1	1	X-Thick	3.4 ± 0.6	142201	<input type="checkbox"/>	
	AlloDerm SELECT™ Regenerative Tissue Matrix  Non-Meshed Rectangle	Per Request	Per Request	X-Thin	0.55 ± 0.25	141400	<input type="checkbox"/>
Per Request		Per Request	Thin	1.0 ± 0.2	141600	<input type="checkbox"/>	
AlloDerm SELECT™ Regenerative Tissue Matrix  Meshed Rectangle	5 x 6	30	X-Thin	0.55 ± 0.25	141430M	<input type="checkbox"/>	
	Per Request	Per Request	X-Thin	0.55 ± 0.25	141400M	<input type="checkbox"/>	
AlloDerm SELECT™ Regenerative Tissue Matrix  Fenestrated Rectangle	4 x 4	16	Thin	1.0 ± 0.2	141616F	<input type="checkbox"/>	
	4 x 8	32	Thin	1.0 ± 0.2	141632F	<input type="checkbox"/>	

SHIPPING REQUIREMENTS: Orders placed by 7:00 PM ET, Monday through Thursday, will ship same day for next-day delivery. AlloDerm SELECT™ Regenerative Tissue Matrix is shipped at ambient temperature via overnight delivery service within 24 hours of order placement.

STORAGE: Allergan products must be stored according to label.

To place an order,
please call Allergan at
1.800.367.5737 or fax
form to 1.800.226.2714.